**POSTURE
TO MAKE YOUR
ENERGY LAST**

**USE
THE RIGHT MUSCLES
WHEN LIFTING**

**POSTURE
REMINDS YOU TO USE
YOUR OWN MUSCLES**

**ONLY WEAR YOUR
POSTURE A FEW HOURS
TO START WITH**

**PERFECT POSTURE
WHILE PERFORMING**

1.

Open the velcro leash and set the regulatory system to desired position. Another person nearby could help you with this after step 2. Close the velcro leash afterwards.

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
Slip in one arm at a time in the loops and position the brace comfortably over the shoulders.

The adjustable strap should sit just below the shoulder blades.

Adjust the velcro strap to increase or decrease the corrective force.

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

BLACKROLL®
THE HEALTHSTYLE COMPANY

POSTURE 2.0

A KIND REMINDER TO STRAIGHTEN UP

BLACKROLL AG
Hauptstrasse 17
CH-8598 Bottighofen

WWW.BLACKROLL.COM

EN Sitting in front of a computer all day can have negative consequences: Your shoulders roll forward and your spine becomes excessively curved in the thoracic region. The BLACKROLL® POSTURE 2.0 prevents this by reminding you to keep your shoulders back and your spine properly aligned. The textile material provides gentle encouragement but not too much support so that your muscles still have to do the work and become stronger in the process.

PUTTING IT ON

- Wear the POSTURE 2.0 over a T-Shirt or any other base layer.
- First adjust the regulatory system on the back to desired size/intensity and close the velcro leash.
- Put the POSTURE 2.0 on like a jacket by slipping the arms through the openings.
- Adjust the size with the velcro straps.

THE RIGHT FIT

- Tighten the POSTURE 2.0 that two fingers can fit between your shoulder and the strap.
- You purposely won't find a sizescale on the POSTURE 2.0. Adjust the POSTURE 2.0 in a way you feel comfortable wearing it. In any case it should not feel uncomfortable.
- When you are standing upright, you should barely feel the POSTURE 2.0.
- When your shoulders are slouching forward, you should feel a slight pull. This is a reminder to stand tall and pull your shoulders backward.

THE RIGHT TIME

- Start by wearing the POSTURE 2.0 for short periods of time.
- Increase the time when your muscles are getting stronger.
- Note that wearing the POSTURE 2.0 should never feel uncomfortable.

WHERE TO WEAR IT

- At home and at the office.
- At the gym to maintain your posture during exercise.
- While running or walking and after physical exercise.

CONTRAINDICATIONS

Please consult your doctor or therapist if you suffer from back or neck pain or other injuries before using the POSTURE 2.0. The POSTURE 2.0 is not suited for children under 12 years of age.

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